

# Biotinylated Recombinant Human IL-2 R beta/CD122 Protein

RPCB0525



## Product Information

**Product SKU:** RPCB0525      **Gene ID:** 3560      **Size:** 100µg  
**Tag:** C-His&Avi      **Reactivity:** Human

## Additional Information

**Expression Host:** -      **Swissprot:** P14784  
**Purity:** > 95% by SDS-PAGE ; > 95% by HPLC

## Protein Information

**Background:** IL-2 R beta is a member of the cytokine receptor superfamily. Human IL-2 R beta cDNA encodes a 551 amino acid residue precursor Type I membrane protein with a 26 residue signal peptide, a 214 residue extracellular region, a 25 residue transmembrane region and a 286 residue cytoplasmic domain. Functional IL-2 receptors can exist in two affinity states on cell surfaces, the high affinity complex consisting of heterotrimers of the alpha, beta, and gamma chains, and the intermediate affinity complex comprising heterodimers of the beta and gamma chains.

**Protein Description:** High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human IL-2 R beta/CD122 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method

**Formulation:** Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.