

## RPCB0555

### Product Information

<b>Product SKU:</b>	RPCB0555	<b>Gene ID:</b>	4680	<b>Size:</b>	100µg
<b>Tag:</b>	C-His&Avi	<b>Reactivity:</b>	Human		

### Additional Information

<b>Expression Host:</b>	-	<b>Swissprot:</b>	P40199
<b>Purity:</b>	> 95% by SDS-PAGE ; > 95% by HPLC		

### Protein Information

**Background:** Carcinoembryonic antigen-related cell adhesion molecule 6 (CEACAM6) belongs to the human carcino-embryonic antigen (CEA) family. Numerous lines of studies have indicated that altered expression of CEACAM6 may have a role in carcinogenesis and development.

**Protein Description:** High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human CEACAM6/CD66c Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method

**Formulation:** Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.