

Recombinant Human Fc gamma RI/FCGR1/CD64 Protein

RPCB0576



Product Information

Product SKU:	RPCB0576	Gene ID:	2209	Size:	100µg
Tag:	C-His&Avi	Reactivity:	Human		

Additional Information

Expression Host:	E. coli	Swissprot:	P12314-1
Purity:	> 95% by SDS-PAGE ; > 95% by HPLC		

Protein Information

Background: The Fc gamma Rs have been divided into three classes based on close relationships in their extracellular domains; these groups are designated Fc gamma¹ RI (also known as CD64), Fc gamma² RII (CD32), and Fc gamma³ RIII (CD16). Each group may be encoded by multiple genes and exist in different isoforms depending on species and cell type.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Fc gamma RI/FCGR1/CD64 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.