

RPCB0610

---

## Product Information

<b>Product SKU:</b> RPCB0610	<b>Gene ID:</b> 3601	<b>Size:</b> 100µg
<b>Tag:</b> C-hFc	<b>Reactivity:</b> Human	

---

## Additional Information

<b>Expression Host:</b> HEK293 cells	<b>Swissprot:</b> Q13261-1
<b>Purity:</b> > 95% by SDS-PAGE ; > 95% by HPLC	

---

## Protein Information

**Background:** Interleukin-15 receptor alpha (IL-15R alpha) is a high affinity IL-15 binding protein that is crucial for mediating IL-15 functions such as memory CD8 T cell proliferation and NK, NK/T cell, and intestinal intraepithelial lymphocyte development.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human IL-15RA/CD215 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method

**Formulation:** Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.