

Biotinylated Recombinant Human IL-12/IL-12A&IL-12B Protein

RPCB0664



Product Information

Product SKU: RPCB0664 **Gene ID:** 3592&3593 **Size:** 100µg
Tag: C-His&Avi **Reactivity:** Human

Additional Information

Expression Host: - **Swissprot:** P29459&P29460
Purity: > 95% by SDS-PAGE ; > 95% by HPLC

Protein Information

Background: Interleukin 12 (IL-12) is a pleiotropic cytokine that plays an essential role in Th1-type immune response against cancer, a condition where cells in a particular part of the body grow and reproduce uncontrollably.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human IL-12/IL-12A&IL-12B Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.