

Biotinylated Recombinant Human CD3E&CD3G Protein

RPCB0679



Product Information

Product SKU:	RPCB0679	Gene ID:	916&917	Size:	100µg
Tag:	C-hFc&Avi	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	P07766(CD3E)&P09693(CD3G)
Purity:	> 95% by SDS-PAGE ; > 95% by HPLC		

Protein Information

Background: T-cell surface glycoprotein CD3 epsilon&CD3 gamma chain, also known as CD3E&CD3G, are single-pass type I membrane proteins. When antigen presenting cells (APCs) activate T-cell receptor (TCR), TCR-mediated signals are transmitted across the cell membrane by the CD3 chains CD3D, CD3E, CD3G and CD3Z. All CD3 chains contain immunoreceptor tyrosine-based activation motifs (ITAMs) in their cytoplasmic domain.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human CD3E&CD3G Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.