

**RPCB0750**

---

## Product Information

<b>Product SKU:</b>	RPCB0750	<b>Gene ID:</b>	64115	<b>Size:</b>	10µg
<b>Tag:</b>	C-hFc&His	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	AAH20568.1
<b>Purity:</b>	> 97% by SDS-PAGE.		

---

## Protein Information

**Background:** V-type immunoglobulin domain-containing suppressor of T-cell activation(VISTA) also known as platelet receptor Gi24 , stress-induced secreted protein-1 (Sisp-1) and B7

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human B7-H5/Gi24/VISTA Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.