

# Recombinant Human Angiopoietin-like 3/ANGPTL3 (17-220) Protein

RPCB0779



## Product Information

**Product SKU:** RPCB0779      **Gene ID:** 27329      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Human

## Additional Information

**Expression Host:** -      **Swissprot:** Q9Y5C1  
**Purity:** > 90% by SDS-PAGE.

## Protein Information

**Background:** Angiopoietin-like protein 3 (ANGPTL3) also known as Angiopoietin-related protein 3, Angiopoietin-5 (ANG-5) is a secreted glycoprotein that is structurally related to the angiopoietins. Human ANGPTL3 contains an N-terminal coiled coil domain and a C-terminal fibrinogen-like domain. The fibrinogen C-terminal domain is sufficient to mediate endothelial cell adhesion. ANGPTL3 is principally expressed in the liver and can exert activities related to both angiogenesis and lipid metabolism. ANGPTL3 directly inhibits lipoprotein lipase (LPL) and endothelial lipase (EL), enzymes responsible for hydrolyzing circulating triglycerides and HDL phospholipids.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Angiopoietin-like 3/ANGPTL3(17-220) Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.