

RPCB0824

Product Information

Product SKU: RPCB0824 **Gene ID:** 53833 **Size:** 10µg
Tag: C-Fc **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q6UXL0
Purity: > 95% by SDS-PAGE.

Protein Information

Background: IL20RB and IL20RA (MIM 605620) form a heterodimeric receptor for interleukin-20 (IL20; MIM 605619) (Blumberg et al., 2001 [PubMed 11163236]).

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IL-20RB Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Supplied as a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.