

**RPCB0833**

---

## Product Information

<b>Product SKU:</b>	RPCB0833	<b>Gene ID:</b>	3441	<b>Size:</b>	10µg
<b>Tag:</b>	NO-Tag	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Expression Host:</b>	E. coli	<b>Swissprot:</b>	P05014
<b>Purity:</b>	> 95% by SDS-PAGE.		

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human IFN-alpha 4 Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.