

Recombinant Human Leptin receptor/LEP-R/CD295 Protein

RPCB0839



Product Information

Product SKU: RPCB0839 **Gene ID:** 3953 **Size:** 20µg
Tag: C-hFc&His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** P48357
Purity: > 95% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Leptin receptor/LEP-R/CD295 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of 50mM MES, 100mM NaCl,1mM DTT, pH6.5.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.