

RPCB0867

Product Information

Product SKU:	RPCB0867	Gene ID:	1489678	Size:	10µg
Tag:	C-His	Reactivity:	SARS-CoV		

Additional Information

Expression Host:	E. coli	Swissprot:	P59595
Purity:	> 90% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant SARS-CoV Nucleocapsid Protein , tested reactivity inE. coliand has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Please contact us for more information.

Formulation: Supplied as a 0.22 µm filtered solution in 20mM Tris, 500mM NaCl, 0.1mM EDTA, 20% glycerol, pH8.0.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.