

# Recombinant Human Catenin beta-1/CTNNB1 Protein

## RPCB0907



### Product Information

**Product SKU:** RPCB0907      **Gene ID:** 1499      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Human

### Additional Information

**Expression Host:** E. coli      **Swissprot:** P35222-1  
**Purity:** > 90% by SDS-PAGE.

### Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Catenin beta-1/CTNNB1 Protein, tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 1.0 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of 50mM Tris, 150mM NaCl, pH 8.0. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.