

RPCB0925

Product Information

Product SKU:	RPCB0925	Gene ID:	1489668	Size:	100µg
Tag:	C-His	Reactivity:	SARS-CoV		

Additional Information

Expression Host:	-	Swissprot:	P59594
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant SARS-CoV Spike RBD Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, 300mM NaCl, pH 7.4. or Supplied as a 0.22 µm filtered solution in PBS, 300mM NaCl, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.