

# Recombinant Mouse P-Selectin/SELP/CD62P Protein

## RPCB0967



### Product Information

<b>Product SKU:</b>	RPCB0967	<b>Gene ID:</b>	20344	<b>Size:</b>	10µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	Mouse		

### Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	Q01102
<b>Purity:</b>	> 90% by SDS-PAGE.		

### Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Mouse P-Selectin/SELP/CD62P Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.