

RPCB1034

Product Information

Product SKU:	RPCB1034	Gene ID:	3549	Size:	10µg
Tag:	NO-tag	Reactivity:	Human		

Additional Information

Expression Host:	E. coli	Swissprot:	Q14623
Purity:	> 97% by SDS-PAGE.		

Protein Information

Background: By Similarity2 Publications Indian hedgehog protein N-product The dually lipidated indian hedgehog protein N-product is a morphogen which is essential for a variety of patterning events during development.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IHH Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.