

RPCB1081

Product Information

Product SKU: RPCB1081 **Gene ID:** 43740568 **Size:** 100µg
Tag: C-mFc **Reactivity:** SARS-CoV-2

Additional Information

Expression Host: HEK293 cells **Swissprot:** P0DTC2
Purity: > 85% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant SARS-CoV-2 Spike S1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4 or Supplied as a 0.22 µm filtered solution in PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.