

Recombinant Human Lysyl hydroxylase 2/PLOD2 Protein

RPCB1106



Product Information

Product SKU: RPCB1106 **Gene ID:** 5352 **Size:** 10µg
Tag: C-His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** O00469-1
Purity: > 95% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Lysyl hydroxylase 2/PLOD2 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Supplied as a 0.22 µm filtered solution in PBS, pH 7.4.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.