

RPCB1124

Product Information

Product SKU: RPCB1124	Gene ID: 57264	Size: 10µg
Tag: N-hFC	Reactivity: Mouse	

Additional Information

Expression Host: -	Swissprot: Q99P87
Purity: -	

Protein Information

Background: Resistin is an adipocytokine, which has been studied for its role in insulin resistance and recently in inflammation. The RETN and CAP1 polymorphisms and gene expression may be potential biomarkers for breast cancer risk. Resistin (RETN), recently found to be relevant to inflammation and inflammatory disorders.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse RETN Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.01 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.