

RPCB1144

Product Information

Product SKU:	RPCB1144	Gene ID:	2919	Size:	10µg
Tag:	NO-tag	Reactivity:	Human		

Additional Information

Expression Host:	Pichia	Swissprot:	P09341-1
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: Has chemotactic activity for neutrophils. May play a role in inflammation and exerts its effects on endothelial cells in an autocrine fashion. In vitro, the processed forms GRO-alpha(4-73), GRO-alpha(5-73) and GRO-alpha(6-73) show a 30-fold higher chemotactic activity.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human CXCL1/GRO-alpha Protein, tested reactivity in Pichia and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.