

RPCB1169

Product Information

Product SKU: RPCB1169 **Gene ID:** 50706 **Size:** 10µg
Tag: C-6His **Reactivity:** Mouse

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q62009-2
Purity: -

Protein Information

Background: Periostin (POSTN), also known as OSF2 (osteoblast specific factor 2), is a heterofunctional secreted extracellular matrix (ECM) protein comprised of four fasciclin domains that promotes cellular adhesion and movement, as well as collagen fibrillogenesis. Postn is frequently overexpressed in various types of human cancers, stimulating metastatic growth by promoting cancer cell survival, invasion and angiogenesis, and can be a useful marker to predict the behavior of cancer.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse Periostin/PN/POSTN Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.