

RPCB1175

Product Information

Product SKU: RPCB1175 **Gene ID:** 17840 **Size:** 10µg
Tag: C-His **Reactivity:** Mouse

Additional Information

Expression Host: HEK293 cells **Swissprot:** P11588
Purity: > 95% by SDS-PAGE.

Protein Information

Background: Binds pheromones that are released from drying urine of males. These pheromones affect the sexual behavior of females.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse MUP-1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Supplied as a 0.22 µm filtered solution in PBS, pH 7.4.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.