

RPCB1177

Product Information

Product SKU:	RPCB1177	Gene ID:	3479	Size:	50µg
Tag:	No tag	Reactivity:	Human		

Additional Information

Expression Host:	Pichia	Swissprot:	P05019-1
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IGF-I Protein , tested reactivity inPichiaand has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.