

Recombinant Human Angiopoietin-4/ANG-4/ANGPT4 (23-503) Protein

RPCB1196



Product Information

Product SKU: RPCB1196 **Gene ID:** 51378 **Size:** 10µg
Tag: C-His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q9Y264
Purity: > 90% by SDS-PAGE.

Protein Information

Background: Human Angiopoietin-4 (Ang-4) (1), alternatively named Ang-3 (2), is a secreted glycoprotein belonging to the angiopoietin family. It has the characteristic structural motifs of angiopoietins including the coiled-coiled domain near the amino-terminus and a fibrinogen-like domain at the C-terminus.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Angiopoietin-4/ANG-4/ANGPT4(23-503) Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.