

Recombinant HumanAmyloid-beta precursor protein/APP (1-687) Protein

RPCB1208



Product Information

Product SKU: RPCB1208 **Gene ID:** 351 **Size:** 10µg
Tag: C-6His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** P05067
Purity: > 97% by SDS-PAGE.

Protein Information

Background: Amyloid β precursor-like protein 2 (APLP2) has been determined to serve an important role in the progression of a number of cancer types. APLP2 expression was significantly associated with disease-specific survival ($P < 0.001$). APLP2 may be used to potentially predict patient prognosis, and to guide clinical diagnosis and treatment in CCRCC.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant HumanAmyloid-beta precursor protein/APP (1-687) Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.