

# Recombinant Mouse Alpha-2-MRAP/LRPAP1 Protein

## RPCB1274



### Product Information

<b>Product SKU:</b>	RPCB1274	<b>Gene ID:</b>	16976	<b>Size:</b>	10µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	Mouse		

### Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	P55302
<b>Purity:</b>	> 95% by SDS-PAGE.		

### Protein Information

**Background:** This protein that interacts with the low density lipoprotein (LDL) receptor-related protein and facilitates its proper folding and localization by preventing the binding of ligands.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Mouse Alpha-2-MRAP/LRPAP1 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** <0.1EU/µg

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.