

RPCB1301

---

**Product Information**

**Product SKU:** RPCB1301      **Gene ID:** 308      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Human

---

**Additional Information**

**Expression Host:** E. coli      **Swissprot:** P08758  
**Purity:** > 97% by SDS-PAGE.

---

**Protein Information**

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Annexin A5/ANXA5 Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 1.0 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.