

**RPCB1361**

---

## Product Information

**Product SKU:** RPCB1361      **Gene ID:** 1803      **Size:** 10µg  
**Tag:** C-His      **Reactivity:** Human

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** P27487  
**Purity:** > 95% by SDS-PAGE.

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human DPP4/CD26 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Supplied as a 0.22 µm filtered solution in PBS, pH 7.4.

**Storage:** Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.