

RPCB1417

---

**Product Information**

**Product SKU:** RPCB1417      **Gene ID:** 102118696      **Size:** 100µg  
**Tag:** C-His      **Reactivity:** Cynomolgus

---

**Additional Information**

**Expression Host:** HEK293 cells      **Swissprot:** G7PG38  
**Purity:** > 95% by Tris-Bis PAGE;> 95% by SEC-HPLC

---

**Protein Information**

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus TNFRSF5/CD40 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.