

RPCB1436

---

## Product Information

**Product SKU:** RPCB1436      **Gene ID:** 102119976      **Size:** 100µg  
**Tag:** C-His&Avi      **Reactivity:** Cynomolgus

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** XP\_005552816.1  
**Purity:** > 95% by Tris-Bis PAGE

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus IL-17A/CTLA-8 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.