

**RPCB1479**

---

## Product Information

<b>Product SKU:</b>	RPCB1479	<b>Gene ID:</b>	43740568	<b>Size:</b>	100µg
<b>Tag:</b>	C-hFc&Avi	<b>Reactivity:</b>	SARS-COV-2		

---

## Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	QHD43416.1
<b>Purity:</b>	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant SARS-COV-2 Spike S1 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.