

**RPCB1504**

---

## Product Information

<b>Product SKU:</b> RPCB1504	<b>Gene ID:</b> 2920	<b>Size:</b> 10µg
<b>Tag:</b> NO-tag	<b>Reactivity:</b> Human	

---

## Additional Information

<b>Expression Host:</b> -	<b>Swissprot:</b> P19875
<b>Purity:</b> > 95% by SDS-PAGE.	

---

## Protein Information

**Background:** Produced by activated monocytes and neutrophils and expressed at sites of inflammation. Hematopoietic progenitor cell proliferation. GRO-beta(5-73) shows a highly enhanced hematopoietic activity.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human CXCL2/MIP-2 Protein, tested reactivity in Pichia and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** <0.1EU/µg

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.