

RPCB1510

Product Information

Product SKU:	RPCB1510	Gene ID:	3592&3593	Size:	10µg
Tag:	C-His	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	P29459&P29460
Purity:	-		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IL-12/IL-12A&IL-12B Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: <0.01EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.