

Recombinant Cynomolgus Siglec-15/CD33L3 Protein

RPCB1519



Product Information

Product SKU: RPCB1519 **Gene ID:** 102126346 **Size:** 10µg
Tag: C-hFc&His **Reactivity:** Cynomolgus

Additional Information

Expression Host: HEK293 cells **Swissprot:** A0A2K5UY47
Purity: > 90% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus Siglec-15/CD33L3 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, 300mM NaCl, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.