

RPCB1551

Product Information

Product SKU:	RPCB1551	Gene ID:	152404	Size:	10µg
Tag:	C-hFc&Avi	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	Q5DX21-1
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: Immunoglobulin superfamily member 11(IGSF11) is expressed on the plasma membrane in the testis and brain. These IGSF proteins undergo final modifications during capacitation and/or the acrosome reaction. IGSF proteins share significant homology with endothelial cell-selective adhesion molecule and coxsackievirus and adenovirus receptor, which mediates cell attachment and homotypic intercellular interactions. In clinical, the IGSF11 has been reported to overexpressed in colorectal cancers and hepatocellular carcinomas, as well as intestinal-type gastric cancers, compared to their corresponding non-cancerous tissues. The IGSF11 has also been found expressed abundantly in the testis and ovary and the IGSF11 can be used as a candidate of cancer-testis antigen.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human VSIG3/IgSF11 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.