

Product Information

Product SKU:	RPCB1694	Gene ID:	8807	Size:	100µg
Tag:	C-hFc	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	O95256-1
Purity:	> 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC		

Protein Information

Background: lactoferrin (LTF), peptidoglycan recognition protein 1 and IL-18 receptor accessory protein (IL18RAP) remained significantly elevated in peripheral monocytes of hypertensive individuals in a separate validation cohort. Interestingly, IL18RAP expression associated with MAP in a cohort of African Americans..

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IL-18 R beta/IL-1 R7/CD218b Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.