

RPCB1730

Product Information

Product SKU:	RPCB1730	Gene ID:	21944	Size:	50µg
Tag:	C-raFc	Reactivity:	Mouse		

Additional Information

Expression Host:	HEK293 Cells	Swissprot:	O54907
Purity:	> 90% by SDS-PAGE.		

Protein Information

Background: TNFSF12 is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. It is a ligand for the FN14/TWEAKR receptor. TNFSF12 has overlapping signaling functions with TNF, but displays a much wider tissue distribution. It can induce apoptosis via multiple pathways of cell death in a cell type-specific manner. It is also found that TNFSF12 promotes proliferation and migration of endothelial cells, and thus acts as a regulator of angiogenesis. TNFSF12 also is a weak inducer of apoptosis in some cell types and mediates NF-kappa-B activation.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Mouse TNFSF12/TWEAK Protein, tested reactivity in HEK293 Cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.