

RPCB1744

Product Information

Product SKU: RPCB1744 **Gene ID:** 3601&3600 **Size:** 100µg
Tag: C-His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** Q13261-1(IL-15RA)&P40933-1(IL-15)
Purity: > 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC

Protein Information

Background: Interleukin-15 receptor alpha (IL-15R alpha) is a high affinity IL-15 binding protein that is crucial for mediating IL-15 functions such as memory CD8 T cell proliferation and NK, NK/T cell, and intestinal intraepithelial lymphocyte development.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IL-15RA&IL-15 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.