

RPCB1746

Product Information

Product SKU:	RPCB1746	Gene ID:	8788	Size:	100µg
Tag:	C-His	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P80370-1
Purity:	> 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC		

Protein Information

Background: paternally inherited genetic defects of DLK1 were identified in four families with nonsyndromic CPP and a metabolic phenotype. DLK1 encodes a transmembrane protein that is important for adipose tissue homeostasis and neurogenesis and is located in the imprinted chromosome 14q32 region associated with Temple syndrome.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human DLK1 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.