

Product Information

Product SKU:	RPCB1758	Gene ID:	-	Size:	100µg
Tag:	C-His	Reactivity:	Cynomolgus		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	XP_005590346.1
Purity:	> 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC		

Protein Information

Background: LILRA4, also known as ILT7 and CD85g, is an approximately 60-70 kDa variably glycosylated transmembrane protein that regulates immune cell activation. Mature human LILRA4 consists of a 423 amino acid (aa) extracellular domain (ECD) with four immunoglobulin-like domains, a 21 aa transmembrane segment, and a 32 aa cytoplasmic domain. LILRA4 function coreceptor to limit the innate immune responses to viral infections; signaling occurs via FCER1G (PubMed:16735691, PubMed:19564354). Down-regulates the production of IFNA1, IFNA2, IFNA4, IFNB1 and TNF by plasmacytoid dendritic cells that have been exposed to influenza virus or cytidine-phosphate-guanosine (CpG) dinucleotides, indicating it functions as negative regulator of TLR7 and TLR9 signaling cascades.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus LILRA4/ILT-7/CD85g Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.