

Recombinant Human Alkaline phosphatase (Intestinal type)/ALPI Protein

RPCB1798



Product Information

Product SKU: RPCB1798 **Gene ID:** 248 **Size:** 100µg
Tag: N-His&Flag **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** P09923
Purity: > 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC

Protein Information

Background: Alkaline phosphatase can be considered "our favorite enzyme" for reasons apparent to those who diagnose and treat metabolic bone diseases or who study skeletal biology. Few might know, however, that alkaline phosphatase likely represents the most frequently assayed enzyme in all of medicine. Elevated activity in the circulation is universally recognized as a marker for skeletal or hepatobiliary disease.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Alkaline phosphatase (Intestinal type)/ALPI Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.