

**RPCB1895**

---

## Product Information

**Product SKU:** RPCB1895      **Gene ID:** 10437      **Size:** 100µg  
**Tag:** C-His      **Reactivity:** Human

---

## Additional Information

**Expression Host:** HEK293 cells      **Swissprot:** P13284  
**Purity:** > 96 % as determined by SDS-PAGE

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human IFI30 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method

**Formulation:** Lyophilized from sterile PBS, pH 7.4

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.