

Recombinant Cynomolgus B7-H1/PD-L1/CD274 Protein

RPCB1910



Product Information

Product SKU:	RPCB1910	Gene ID:	102145573	Size:	100µg
Tag:	C-hFc	Reactivity:	Cynomolgus		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	G7PSE7
Purity:	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus B7-H1/PD-L1/CD274 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.