

RPCB1959

Product Information

Product SKU:	RPCB1959	Gene ID:	10397	Size:	100µg
Tag:	C-His	Reactivity:	Human		

Additional Information

Expression Host:	E.coli	Swissprot:	Q92597
Purity:	> 85% by SDS-PAGE.		

Protein Information

Background: NDRG1 is a cytoplasmic protein involved in stress responses, hormone responses, cell growth, and differentiation. NDRG1 is necessary for p53-mediated caspase activation and apoptosis. Mutations in the NDRG1 gene are a cause of Charcot-Marie-Tooth disease type 4D, and expression of this gene may be a prognostic indicator for several types of cancer.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human NDRG1 Protein, tested reactivity in E.coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: Please contact us for more information.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.