

Recombinant Human Follistatin-like protein 1/FSTL1 Protein

RPCB2012



Product Information

Product SKU:	RPCB2012	Gene ID:	11167	Size:	10µg
Tag:	C-6*His	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q12841-1
Purity:	PBS		

Protein Information

Background: Follistatin-like 1 (FSTL1), a secreted glycoprotein, has been shown to participate in regulating developmental processes and to be involved in states of disease and injury. Recent findings on FSTL1 in both acute coronary syndrome and heart failure emphasize its potential as a target for cardiac regenerative therapy.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Follistatin-like protein 1/FSTL1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.01 EU/µg

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.