

RPCB2025

Product Information

Product SKU:	RPCB2025	Gene ID:	141	Size:	100µg
Tag:	N-His	Reactivity:	Human		

Additional Information

Expression Host:	E.coli	Swissprot:	P54922
Purity:	> 90% by SDS-PAGE.		

Protein Information

Background: ADPRH belongs to the ADP-Ribosylglycohydrolase family. ADPRH specifically acts as an arginine mono-ADP-ribosylhydrolase by mediating the removal of mono-ADP-ribose attached to arginine residues on proteins. These reactions are related to cell signaling and the control of many cell processes, such as DNA repair and cell apoptosis.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human ADPRH Protein, tested reactivity in E.coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: Please contact us for more information.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.