

Product Information

Product SKU: RPCB2127 **Gene ID:** 3001 **Size:** 10µg
Tag: C-His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** P12544
Purity: >95% as determined by SDS-PAGE

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Granzyme A/CTLA-3/GZMA Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.