

RPCB2132

Product Information

Product SKU: RPCB2132 **Gene ID:** 998 **Size:** 100µg
Tag: N-GST **Reactivity:** Human

Additional Information

Expression Host: E.coli **Swissprot:** P60953-2
Purity: > 85% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human CDC42 Protein , tested reactivity in E.coli and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: -

Formulation: Lyophilized from sterile 20mM Tris, 0.15M NaCl, 0.5mM GSH, pH 8.0. Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.