

Protein

RPCB2178

Product Information

Product SKU: RPCB2178 **Gene ID:** - **Size:** 100µg
Tag: C-His&Avi **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** A0A140T913(HLA-A*02:01)&P61769(B2M)&SLLMWITQV
Purity: > 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC

Protein Information

Background: NY-ESO-1 or New York esophageal squamous cell carcinoma 1 is a well-known cancer-testis antigen (CTAs) with re-expression in numerous cancer types. Its ability to elicit spontaneous humoral and cellular immune responses, together with its restricted expression pattern, have rendered it a good candidate target for cancer immunotherapy.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human NY-ESO-1 (HLA-A*02:01) Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.