

RPES8079

Product Information

Product SKU:	RPES8079	Expression Host:	E.coli	Size:	20µg
Tag:	N-Trx	Reactivity:	Human	Accession:	P20702

Additional Information

Calculated MW:	56.4 kDa	Observed MW:	56 kDa
Sequence:	Asn21-Ala352		

Protein Information

Background: -

Synonyms: SLEB, Leukocyte adhesion receptor p150+, ITGAX, CD11C, SLEB6, CD11 antigen-like family member C, Integrin alpha-X, Leukocyte adhesion receptor p150+95, p150-95, CR4, integrin α , alpha-X-X

Endotoxin: < 10 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.